**DENOMINAZIONE / PRODUCT
NAME / NOM DU PRODUIT /
PRODUKTBEZEICHNUNG /
DENOMINACIÓN / NOME DO
PRODUTO**

Eurosteril®

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Eurosteril® sterilization pouches and rolls are used for the sterilization of medical-surgical devices and of materials for medical use which are especially suitable for steam or EtO GAS. They are made by 2 layers: one made of high resistance medical paper and the other made of blue polyester/polypropylene. They are available as prefolded and self-adhesive rolls and pouches of various sizes in order to satisfy different needs of packaging. Sealing temperature between 170 °C and 200 °C The seal is made by two parallel strings which are separated by a groove for rolls and sterilization pouches with a width of 5 cm. The seal is made by three parallel strings which are separated by two grooves for rolls and sterilization pouches with other dimensions. On the paper layer two chemical indicators of process are printed: one for steam and the other for ethylene-oxyde, both indicating that the sterilization process has been carried out.

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Eurosteril® pouches and rolls are intended to contain medical devices, resist the sterilization process allowing the sterilization of the contents, allow the preservation of the contents in sterile conditions.

Eurosteril® pouches and reels are compatible with the following sterilization treatments:

- STEAM
- EtO (Ethylene Oxide)

Pouches and reels are equipped with Type 1 process indicators (EN ISO 11140-1) for Ethylene Oxide and Vapor, printed with water-based ink. These provide an indication for the user as they identify whether the barrier system, together with its contents, have already undergone one of the sterilization treatments indicated above. Indicators follow the criterion:

- Steam sterilization = STEAM indicator from pink to brown
- Sterilization with ethylene oxide = EO indicator from blue to yellow / brown

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The sealing on the 4th side for pouches, and on the 3rd and 4th side for reels, can be carried out with all thermosealing machines suitable for the construction materials of Eurosteril® products. It is advisable to refer to the International Standard EN 868-5 for the correct setting parameters (temperature, pressure, speed) in order to guarantee the desired minimum resistance according to the sterilization method employed.

The melting temperature of polypropylene is 175°C.

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Medical devices must be disposed of in accordance with medical procedures adopted by health care institutions and organizations and local, state or federal laws and regulations.

Sterile Barrier System devices of the Eurosteril® series are comparable to urban waste and therefore they must be disposed of as urban waste (mixed, non-recyclable).

The PET plastic film has been processed in such a way as to obtain a stratification that makes it non-recyclable. It must therefore be disposed of as non-recyclable waste.

If the operator decouples film from paper, the paper can be recycled in separate collection.

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