Medi-Test Mikroalbumin

de

In-vitro-Diagnostikum

en

Medi-Test Microalbumin

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Test strips for rapid detection of albumin and creatinine in urine

Intended Use

Medi-Test Microalbumin is a test strip for semi-quantitative determination of albumin and creatinine in urine, as well as for determination of the albumin/creatinine ratio. As a screening test, it serves for detection of microalbuminuria. Positive results indicate incipient kidney disease or complications of cardiovascular disease with increased mortality.¹ Patients at increased risk for developing kidney disease are primarily patients with diabetes and /or hypertension. The test strips are suitable for visual evaluation based on the color scale and table or for automated analysis on the URYXXON® Relax device.

Application only by professional staff.

Instructions for Use

Use fresh, uncentrifuged urine (spot urine) and mix the urine sample well. The urine should be examined within no more than 2 hours of collection. Withdraw a test strip and seal the container again. Dip the test strip into the urine for about 1 second. Wipe off the edge on a paper towel to remove excess urine. In case of visual evaluation, compare the reaction colors after 60 seconds to the color scale and assign to the appropriate value. Subsequently the albumin /creatinine ratio is determined using the interpretation table. Color changes that occur after elapse of the reaction time may not be considered. Additionally, the test strip can be evaluated using the URYXXON® Relax device. To this end, insert the test strip into the device as described in the operating instructions. The test fields are evaluated by reflection photometry, and the results are printed onto the report form. The output values correspond to the visual concentration range, and the albumin / creatinine ratio is assessed in analogy to the interpretation table. In the evaluation it should be noted that in some cases the URYXXON® Relax results may differ from those of the visual evaluation, because the human eye may perceive colors differently than the URYXXON® Relax does.

Principle

Albumin: The test is based on the principle of the “protein error” of indicators, i.e. at a constantly buffered pH, albumin reacts with a tetrabromophenolsulphonephthaleine derivative from yellow-green to green-blue.

Creatinine: The detection is based on the reaction of creatinine with dinitrobenzoic acid. The resulting coloration runs from yellow-brown to blue-black.

Assessment and sources of error

Albumin /creatinine ratio and interpretation table: ²

The test detects values from 10 mg/L albumin and 10 mg/dL creatinine.

Measurement ranges

10 – 150 mg/L albumin
10 – 300 mg/dL creatinine

Based on the albumin and creatinine values determined, the urine sample is assessed in accordance with the value ranges of the interpretation table, either visually (see container label) or automated using the URYXXON® Relax. This method allows use of arbitrary urine samples (please refer to section “Interferences“!) independent of the concentration of the urine. Thus no 24-hour urine collection is required. Using the table, the results are assigned to the urinary classifications of “Normal”, “Abnormal” or “High Abnormal”.

Albumin mg/L	Creatinine mg/dL				
	10	50	100	200	300
10	*				
30					Normal
80					
150					

* Please repeat with a new sample, since it is too strongly diluted.

The albumin / creatinine ratios classified in the table are based on the following ranges of values (mg albumin / g creatinine) ²

Normal: < 30 mg/g

Abnormal: 30 – 299 mg/g (microalbuminuria)

High abnormal: ≥ 300 mg/g (macroalbuminuria or proteinuria)

Reference values

Albumin: In normalbuminuria, albumin concentrations in urine are below 20 mg/L.³

Creatinine: Creatinine is normally present in urine in a concentration range of 10 to 300 mg/dL.

Assessment example

If a urine sample contains 30 mg/L albumin and 200 mg/dL creatinine, the sample is to be classified as NORMAL according to the table. If the test yields a result of 80 mg/L albumin and 100 mg/dL creatinine, the urine sample is to be assessed as ABNORMAL (microalbuminuria).

Interferences

Albumin excretion into urine may be temporarily increased after physical exertion, after urinary tract infections and during febrile illness. Any contamination of the urine samples, vessels or working surfaces by cleaning agents, disinfectants, or urine preservatives (e.g. quaternary ammonium compounds) can result in distorted albumin and creatinine values (boric acid up to 1 g/L has no impact!). Intensely colored urinary components such as haemoglobin (≥ 5 mg/dL) or bilirubin (≥ 2 mg/dL) also interfere with the determination, since the natural color of such urine may obscure the reaction colors. Thus already intensely dark-colored urine is not suitable for the test procedure. False-positive results may occur with highly alkaline urine (pH > 9). Vitamin C does not interfere with the determination.

Reacting substances

Amount per cm² for impregnation

Albumin: Tetrabromophenolsulphonephthaleine derivative 23 µg
Creatinine: Dinitrobenzoic acid 2586 µg

Quality Control

The test strips are checked using commercially available control solutions. The controls should be performed once a day, after opening of a new container of test strips, whenever starting use of a new lot, and every 30 days for verification of the storage conditions. The quality control data must match the values of the acceptance checks. If these are not met, no patient samples may be tested until the problem is resolved. The quality control tests must be repeated until acceptable results are obtained.

Performance Data

For determination of the performance data of Medi-Test Microalbumin, urine samples from urological practice were used which had been obtained as part of daily routine urine control. The measured albumin and creatinine values were compared with the results of the respective reference method (immunoturbidimetric assay for albumin and colorimetric Jaffé method for creatinine).

	negative percent agreement	positive percent agreement	overall percent agreement
Visually (N=102)			
Albumin:	91.9 %	92.5 %	92.2 %
Albumin /creatinine ratio:	92.5 %	88.6 %	91.2 %
URYXXON® Relax (N=612)			
Albumin:	83.6 %	85.0 %	84.2 %
Albumin /creatinine ratio:	84.8 %	89.5 %	86.4 %

Here a negative assessment corresponds to the result “Normal“, a positive assessment to the result “Abnormal“ or “High Abnormal“.

Notes

Since the albumin content in the urine varies, it is recommended to measure several urine samples within 3–6 months. Basically, individual test strip results can allow definitive diagnosis and targeted therapy only in combination with other medical findings. All positive findings should be confirmed by a quantitative method for the determination of albumin. The effect of drugs or their metabolites on the test is not known in all cases. In case of doubt, it is therefore recommended to repeat the test after discontinuation of the medication. For urine collection use only well-washed and clean vessels. Boric acid as a urine preservative does not disturb the test (see above). Always withdraw only the required number of test strips. Immediately reseal package tightly after withdrawal of strips. Do not touch reaction areas!

Storage

Avoid exposing the strips to sunlight and moisture. Keep container in a cool, dry place (storage temperature not above 30 °C /86 °F). If stored properly, the test strips can be used until the printed expiration date.

Important: Before use, make sure that the test strips and container are not damaged, and the expiry date and /or the max. storage temperature have not been exceeded. In these cases, container and test strips must be discarded.

The plug of the test strip container contains a non-toxic desiccant. If it should be accidentally swallowed, drink plenty of water.

For an explanation of the symbols, please see the end of the package insert.

Disposal: Dispose of the used test strips in accordance with the applicable safety regulations.

Commercial form: Package of 24 test strips

Date of revision: 06/2022

References

¹ Apostolovic, S.; Stanojevic, D.; Djordjevic, V.; Jankovic Tomasevic, R.; Salinger Martinovic, S.; Djordjevic Radojkovic, D.; Koracevic, G.; Stojanovic, I.; Pavlovic, M. Clin. Lab. 57; 229-235; 2011

² Position Statement: Nephropathy in Diabetes. *Diabetes Care*, 27, S79-S83 (Supplement 1), 2004

³ Lambers Heerspink, H. J.; Brantsma, A. H.; de Zeeuw, D.; Bakker; S. J. L.; de Jong; P. E.; Gansevoort; R. T. *Am. J. Epidemiol.* 168; 897-905; 2008

	Hersteller / Manufacturer / Fabricante / Fabricant
	Artikelnummer / Item number / Referencia / Référence produit
	Gebrauchsanweisung beachten / Please read instructions for use! / Obsérvense las instrucciones de uso. / Respecter les instructions d'utilisation
	Temperaturbegrenzung / Permitted storage temperature range / Límites de temperatura / Limites de température
	Nicht wiederverwenden / Do not reuse / Producto de un solo uso / Ne pas réutiliser

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany
Tel.: +49 24 21 969-0 · info@mn-net.com · **www.mn-net.com**

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