TMS Uri STIX 10

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In-vitro-Diagnostikum

EN

Test strips for rapid determination of blood, urobilinogen, bilirubin, protein, nitrite, ketones, ascorbic acid, glucose, pH-value and leukocytes in urine

Use

Screening test for detection of diabetes, metabolic abnormalities, liver diseases, biliary and hepatic obstructions, hemolytic diseases and diseases of kidney and urinary tract.

Only for use by qualified personnel.

Instructions for use

Dip the test strip for approximately 1 second into the fresh urine. Draw it across the rim of the container to remove excess urine. After 30 to 60 seconds (leukocyte test field after 60–120 seconds) compare the test strip with the color scale. The best time for comparison is after 30 seconds. Color changes that take place after more than 2 minutes are of no significance. When tested, the urine should not be older than 2 hours.

Principle

Blood: The detection is based on the pseudoperoxidative activity of hemoglobin and myoglobin, which catalyze the oxidation of an indicator by an organic hydroperoxide producing a green color.

Urobilinogen: The test paper contains a stable diazonium salt, forming a reddish azo compound with urobilinogen.

Bilirubin: A red azo compound is obtained in the presence of acid by coupling of bilirubin with a diazonium salt.

Protein: The test is based on the “protein error” principle of indicators. The test zone is buffered to a constant pH value and changes color from yellow to greenish blue in the presence of albumin. Other proteins are indicated with less sensitivity.

Nitrite: Microorganisms, which are able to reduce nitrate to nitrite, are indicated indirectly by this test. The principle of Griess reagent is the basis of this test. The test paper contains an amine and a coupling component. A red colored azo compound is formed by diazotisation and subsequent coupling.

Ketones: The test is based on the principle of Legal’s test. Acetoacetic acid and acetone form with sodium nitroprusside in alkaline medium, a violet colored complex.

Ascorbic acid: The detection is based on the decoloration of Tillmans reagent. In the presence of ascorbic acid, a color change takes place from blue to red.

Glucose: The detection is based on the glucoseoxidase-peroxidase-chromogen-reaction. Apart from glucose, no other compound in urine is known to give a positive reaction.

pH: The test paper contains indicators which clearly change color between pH 5 and pH 9 (from orange to green to turquoise).

Leukocytes: The test is based on the esterase activity of granulocytes. This enzyme splits carboxylic acid ester. The alcohol constituent released reacts with a diazo salt, producing a violet color.

Evaluation – Sources of Error

Blood: The minimum sensitivity of the test strip is 5 erythrocytes/µL urine, corresponding to approx. 0.015 mg hemoglobin/dL urine. Intact erythrocytes are indicated by flecky discolorations of the test field. The color fields correspond to the following values:

0 (negative), ca. 5–10, ca. 50, ca. 250 Ery/µL resp.

hemoglobin concentration out of ca. 10, ca. 50, ca. 250 Ery/µL

Normal concentrations of ascorbic acid (< 40 mg/dL) do not influence the test results. Falsely positive reactions can be produced by a residue of peroxide containing cleansing agents.

Urobilinogen: In dependence upon the urine color, 0,5 to 1 mg urobilinogen/dL urine are indicated. 1 mg/dL is considered to be the normal excretion rate. Higher values are pathological. A complete absence of urobilinogen in the urine, which is likewise pathological, cannot be indicated by the strips. The color fields correspond to the following urobilinogen concentrations:

norm. (normal), 2, 4, 8, 12 mg/dL or

norm. (normal), 35, 70, 140, 200 µmol/L

The test will be inhibited by higher concentrations of formaldehyde. Exposure of the urine to light for a longer period of time may lead to lowered or falsely negative results. Too high or falsely positive results can be caused by the presence of diagnostic or therapeutic dyes in the urine. Larger amounts of bilirubin produce a yellow coloration.

Bilirubin: The minimum sensitivity of the test strip is 0,5 to 0,75 mg bilirubin/dL urine. The color fields correspond to the following values:

0 (negative), 1(+), 2(++), 4(+++) mg/dL or

0 (negative), 17(+), 35(++), 70(+++) µmol/L

Some urine contents can produce a yellow coloration of the test strip. Ascorbic acid and nitrite in higher concentrations inhibit the test. Exposure of the urine to light for a longer period of time may lead to lowered or falsely negative results. Too high or falsely positive results can be caused by the presence of diagnostic or therapeutic dyes in the urine.

Protein: The minimum sensitivity of the test strip is 10 mg protein/dL urine. The color fields correspond to the following ranges of albumin concentrations:

negative, 30, 100, 500 mg/dL or

negative, 0,3, 1,0, 5,0 g/L

Falsely positive results are possible in alkaline urine samples (pH > 9), after infusions with polyvinylpyrrolidone (blood substitute), after intake of medicaments containing quinine and also by disinfectant residues in the urine sampling vessel. The protein coloration may be masked by the presence of medical dyes (e.g. methylene blue) or beetroot pigments.

Nitrite: The test detects concentrations from 0,05 mg nitrite/dL urine. Every pink color indicates a bacterial infection of the urinary tract. The color intensity depends only on the nitrite concentration, but does not provide information about the extent of the infection. A negative result does not preclude an infection of the urinary tract, if bacteria which cannot produce nitrite are present. Falsely negative results can be produced by high doses of ascorbic acid, by antibiotics therapy and by very low nitrate concentrations in urine as the result of low nitrate diet or strong dilution (diuresis). Falsely positive results can be caused by the presence of diagnostic or therapeutic dyes in the urine.

Ketones: The test is more sensitive to acetoacetic acid than to acetone. Values of 5 mg/dL acetoacetic acid or 50 mg/dL acetone are indicated. The color fields correspond to the following acetoacetic acid values:

0 (negative), 25(+), 100(++), 300(+++) mg/dL or

0 (negative), 2,5(+), 10(++), 30(+++) mmol/L

Phenylketones in higher concentrations interfere with the test, and will produce variable colors.

β-Hydroxybutyric acid is not detected. Phthalein compounds interfere by producing a red coloration.

Ascorbic acid: The color fields correspond to the following values:

0 (negative), 10(+), 20(++) mg/dL or

0 (negative), 0,6(+), 1,1(++), mmol/L

Only for information!

Glucose: Pathological glucose concentrations are indicated by a color change from green to bluish green. Yellow or greenish test fields should be considered negative or normal. The color fields correspond to the following ranges of glucose concentrations:

neg. (yellow), neg. or normal (greenish), 50, 150, 500, ≥ 1000 mg/dL or

neg. (yellow), neg. or normal (greenish), 2,8, 8,3, 27,8, ≥ 55,5 mmol/L

The influence of ascorbic acid (vitamin C) has been largely eliminated. An inhibitory effect is produced by gentisic acid. Falsely positive reactions can be produced by a residue of peroxide containing cleansing agents.

pH: The pH value of fresh urine of healthy people varies between pH 5 and pH 6. The color scale gives a clear distinction of pH value between pH 5 and pH 9.

Leukocytes: The test records values starting from approx. 10 leukocytes/µL urine. Changes in color that can not be assigned to the negative reference field and faint violet colors after 120 seconds must be evaluated as positive. The color reference fields correspond to the following leukocyte concentrations:

negative (normal), 25, 75, 500 leukocytes/µL

A weakened reaction can be expected in the case of proteinuria of over 500 mg/dL and a glucose concentration of over 2 g/dL as well as in the case of patients taking preparations containing cephalexin and gentamycin. Bacteria, trichomonads and erythrocytes do not re-act with this test. Formaldehyde (as a preservative) can result in a false positive reaction. Boric acid used as preservative decreases the sensitivity of the reaction. Excretion of bilirubin, nitrofruantoin or other stronglycolored compounds may disguise the color of the reaction. Tests with female patients have shown that vaginal discharge can cause a false positive reaction.

Quality Control in professional use

The performance of the test strips should be confirmed by use of positive and negative control solutions. Positive and negative controls should be analyzed once a day, whenever a new bottle of strips is opened, whenever a new lot of strips is started, and every 30 days to check storage conditions. Each laboratory should establish its own goals for adequate standards of performance, and should question handling and testing procedures if these standards are not met.

Reactive ingredients

(Quantity resp. activity/cm2 at time of impregnation)

Blut :	Nitrite :	Glucose :	
tetramethylbenzidine	31 µg sulfanilic acid	95 µg glucose oxidase	7 U
cumene hydroperoxide	315 µg quinoline derivative	37 µg peroxidase	1 U
Urobilinogen :	Ketones :	tetramethylbenzidine	96 µg
diazonium salt	75 µg sodium nitroprusside	180 µg pH :	
Bilirubin :	Ascorbic acid :	methyl red	3 µg
diazonium salt	29 µg 2,6-dichlorphenolindophenol	7 µg bromothymol blue	10 µg
Protein :		Leukozytes :	
tetrabromophenol blue	10 µg	carboxylic acid ester	16 µg
		diazonium salt	14 µg

Directions

In any case, in order to establish a final diagnosis and prescribe an appropriate therapy, the results obtained with test strips should be verified with other medical results.

The effect of medicament or their metabolic products on the test is not known in all cases. In case of doubt, it is recommended not to take the medicament and then repeat the test. Only use well washed and clean vessels for urine collection. The presence of usual urine preserva-tives will not affect the test results.

Remove only as many test strips as are required, and reseal the container immediately after use. Do not touch the test pads. Avoid exposing the strips to sunlight and moisture. Store the container below + 30 °C in a dry place. The test strips are stable, when stored properly up to the date of expiry indicated.

The caps contain a nonpoisonous and harmless desiccant. In case this desiccant is swallowed accidentally, then drink plenty of water.

Explanation of symbols can be found in the package insert.

Disposal: Please dispose all used test strips in accordance with your local laws and regulations.

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