NADAL[®] Hb/Hp Complex (test cassette)

REF 272011

DE Gebrauchsanweisung	2	PL Sposób użycia	17
EN Instruction for use	5	Symbols	20
FR Instructions d'utilisation	8	Our Teams	20
ES Instrucciones de uso	11		
IT Istruzioni per l'uso	14		

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Friedenstrasse 32
93053 Regensburg
Germany

Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

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14. Referenzen

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; *Archive of Internal Medicine*; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; *British Medical Journal*; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? *Cancer Supplement*; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; *Clinical Biochemistry* (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; *Mobility and Mortality Weekly Report*; (1995) 45 (5) : 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; *Gastroenterology*; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: *Hepatogastroenterology*; 2000 Mar-Apr; 47(32):396-9.

Rev. 1.2, 2013-11-06 OM

1. Intended Use

The NADAL® Hb/Hp Complex is a visual, immunochromatographic dual rapid test for the qualitative detection of human hemoglobin and hemoglobin-haptoglobin complexes in human stool samples.

2. Clinical Significances

The NADAL® Hb/Hp Complex is intended for professional *in vitro* diagnostic use. The test is designed to provide evidence of bowel cancer and/or adenomas. Bowel cancer is one of the most frequently diagnosed types of cancer and one of the most frequent causes of death from cancer. Bowel cancer can be detected early by means of tests for hidden blood in stools and the mortality rate can therefore be reduced. The hemoglobin/haptoglobin (Hb/Hp) complex plays an important role in the retrieval of hemoglobin from lysed erythrocytes. It is relatively stable against acid and proteolytic degradation and can also be detected if it is set free in the anterior intestinal region. Blood from intestinal polyps and higher carcinomas is also detectable. The NADAL® Hb/Hp Complex shows a high sensitivity for the detection of colorectal adenomas and carcinomas. Detection of hemoglobin and hemoglobin-haptoglobin complex is therefore an innovative and extremely sensitive method for the early diagnosis of colon carcinomas and larger adenomas.

3. Principle of the Test

The immunological NADAL® Hb/Hp Complex is based on an immunochromatographic method, in which both Hb and Hb/Hp complexes are detected through specific antibody reactions.

For Hb:

Hemoglobin contained in stool sample reacts with specific gold-conjugated monoclonal antibodies. This complex diffuses along the membrane and binds to the anti-hemoglobin antibody which is immobilised on the test line (T).

For Hb/Hp:

An antibody attaches to hemoglobin and another antibody attaches to haptoglobin. These Hb/Hp complexes are measured and there is no cross-reactivity with hemoglobin. In case of a positive result, hemoglobin and haptoglobin from the stool sample bound to gold-conjugated antibodies bind to immobilised antibodies on the test line (T) and a red/pink colouration becomes visible. In case of a negative result, there are no Hb/Hp complexes that can bind to the antibodies of the test line (T) and no colouration of the test line (T) is formed. A red/pink control line (C) implies that the test has been performed correctly. The test is valid.

4. Reagents and Materials Supplied

- 25 individually packed test cassettes
- 25 sample-taking tubes containing extraction solution
- 25 sets of instructions on taking stool samples
- 25 stool containers
- 1 package insert

5. Additional Required Materials

- 1 stopwatch

6. Storage & Stability

The test can be stored at +2° to +30°C in a sealed bag until the expiration date.

7. Warnings and Precautions

- Only for professional *in vitro* diagnostic use.
- Read the entire package insert carefully prior to testing.
- Do not use the test after the expiration date.
- Do not use the test if the foil pouch is damaged. Do not reuse the test.
- Use a new stool container for each specimen.
- It is recommended to use the enclosed stool container in order to avoid mixing with blood from urine samples or interfering components from the toilet water.
- Do not eat, drink or smoke in any area where specimens and kits are handled. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow standard procedures for the proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Only open the test packaging when the test is actually going to be performed.
- The extraction solution contains a small amount of sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of extraction solution or extracted samples, always flush with copious quantities of water to prevent azide buildup.
- Do not interchange or mix reagents from different lots.
- Humidity and temperature can adversely affect results.
- Used testing materials should be discarded according to local regulations.
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled by observing usual safety precautions (e.g. do not ingest or inhale).
- Patients should follow the instructions for obtaining a sample exactly and should not take any samples during menstruation or 3 days before or after, in case of bleeding hemorrhoids, bleeding caused by constipation or rectal medicine administration.
- Alcohol, aspirin and other medications taken in excess may cause gastrointestinal irritation resulting in occult bleeding. Such substances should be discontinued at least 48 hours prior to testing.
- No dietary restrictions are necessary before testing.

8. Specimen Collection and Preparation

Note: Please give patients the following information.

During sampling, it is important that the sample-taking tube (spiral rod) is inserted **three times** in a row into **three different locations of the stool sample**.

The sample-taking tube should NOT be placed back into the tube between the three collections. No fluid should leak from the tube. **The nature of the test sample and the compliance with instructions will affect the test result.**

Bring the sample-taking tube to room temperature before use. Collect the stool using the stool container enclosed. Hold the sample-taking tube vertically and open the sealing cap. Remove the spiral rod. Place the spiral rod at three different locations in the stool sample. Place the spiral rod with the stool sample into the tube and seal well. Label the well-sealed sample-taking tube with the name of the patient.

Perform testing immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Specimens may be stored at 2-8°C for up to 72 hours.

If specimens are to be shipped, pack them in compliance with all applicable regulations for transportation of etiological agents.

If the sample is taken at the patient's home, the sample should be protected from light and stored in a cool place.

9. Procedure of the Test

1. Bring the test cassette and the patient's sample to room temperature (+20°C to +30°C) before testing.
2. Take the required number of test cassettes from the foil packaging only immediately before performing the test. Mark the cassette with the name of the patient or with any form of identification.
3. Shake the stool sample tube carefully to ensure that the stool sample mixes properly with the extraction solution.
 
4. Unscrew the protective cap. Using a paper towel break the seal of the sample-taking tube with a rotary motion. Hold the sample-taking tube upright and add 3 drops of the solution mixture to each round sample window (S) of the cassette. Avoid trapping air bubbles in the specimen wells (S) and do not add any solution to the result area.
 
5. Read the results after 5 minutes. Strong positive results may be evaluated even sooner. Evaluate the result in a maximum of 15 minutes after dropping the solution mixture into the round sample windows.
 

10. Interpretation of the Results

Test result is evaluated as "positive" if in addition to the control lines at least one of the two test lines develops.


Positive for Hb:

One red/pink line develops in each control line region (C) and one red/pink line develops in the test line region (T) for Hb.

Positive for Hb/Hp:

One red/pink line develops in each control line region (C) and one red/pink line develops in the test line region (T) for Hb/Hp.

Positive for Hb and Hb/Hp:

In addition to the control lines (C), a red/pink line develops in each test line region for Hb and Hb/Hp.

Note: If you are testing strongly positive samples, the intensity of the control line may be reduced. It is not recommended to compare the intensity of the lines.

Negative:

One red/pink line develops in each control line region (C). No line develops in each test line region (T).

Invalid:

The absence of the control lines (C) indicates a failure in the test procedure. The sample should be retested.

Note: Insufficient specimen volume, incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

11. Quality Control

The test contains a procedural control. A red/pink line that develops in each control line region (C) indicates that the test has been performed properly.

The background can become slightly yellowish during testing, depending on the colour of the stool sample. This is acceptable as long as the evaluation of the test results is not adversely affected.

12. Limitations

The test is used for the qualitative detection of human hemoglobin and hemoglobin-haptoglobin complexes in stool samples.

A positive test result will be obtained if human hemoglobin or hemoglobin-haptoglobin complexes are present in the sample. In addition to intestinal bleeding, blood in stools may also have other causes, such as hemorrhoids or a mixture of the stool sample with urine containing blood.

Not all intestinal bleedings are the result of benign or malignant polyps. The information obtained from this test should be examined in conjunction with other clinical test methods.

Carcinomas and polyps can lead to bleeding in different extents and intermittently. It is advisable to check various

stool samples by means of immunological verification procedures.

Urine and excessive dilution of samples with water from the toilet bowl can lead to false test results. It is therefore recommended to use the provided stool container.

Stool samples should not be taken during menstruation or 3 days before or after, in case of bleeding hemorrhoids bleeding caused by constipation or rectal medicine administration. This could lead to false positive results.

Blood-thinning medications, such as for example, acetylsalicylic acid or coumarin as well as iron preparations can lead to bleeding unrelated to tumors.

Unlike test methods that detect only free hemoglobin, the NADAL® Hb/Hp Complex is very sensitive in detecting bleeding in the anterior intestinal passage.

13. Performance Characteristics

Analytical sensitivity

The detection limit of the NADAL® Hb/Hp Complex is 25 ng of free or complex hemoglobin per milliliter of buffer solution.

Prozone effect:

The NADAL® Hb/Hp Complex works reliably at extremely high hemoglobin concentrations (>500,000 ng/mL).

Specificity

The NADAL® Hb/Hp Complex is specific for human hemoglobin or hemoglobin-haptoglobin complex and has no cross-reactivity with hemoglobin from cattle, pigs, rabbits, horses or sheep at a concentration of up to 0.5 mg/mL in the extraction buffer. The NADAL® Hb/Hp Complex also shows no cross-reactivity with Bilirubin, L-ascorbate (vitamin C), and horseradish peroxidase.

14. References

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

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14. Bibliographie

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; *Archive of Internal Medicine*; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; *British Medical Journal*; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, *Why Screen? Cancer Supplement*; (1994) 74(7): 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; *Clinical Biochemistry* (1995) 28(1): 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3): 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; *Morbidity and Mortality Weekly Report*; (1995) 45 (5): 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; *Gastroenterology*; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: *Hepatogastroenterology*; 2000 Mar-Apr; 47(32):396-9.

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14. Referencias

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094

3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

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14. Bibliografia

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; *Archive of Internal Medicine*; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; *British Medical Journal*; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, *Why Screen? Cancer Supplement*; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; *Clinical Biochemistry* (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; *Morbidity and Mortality Weekly Report*; (1995) 45 (5) : 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; *Gastroenterology*; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: *Hepatogastroenterology*; 2000 Mar-Apr; 47(32):396-9.

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nal von minden GmbH