NADAL[®] RSV Test (test cassette)

REF 491905N-20

de	Gebrauchsanweisung	2	cs	Návod k použití	30
en	Instructions for use	6	fi	Käyttöohje	34
fr	Instructions d'utilisation	10	sv	Användarinstruktioner	38
es	Instrucciones de uso	14	no	Bruksanvisning	42
it	Istruzioni per l'uso	18		Symbols	46
pl	Sposób użycia	22		Our Teams	48
pt	Instruções de Utilização	26			

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15. Referenzen

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
3. Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
4. Hall, C,B,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
5. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
6. Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

Rev.0, 2022-05-25 OM/KRM

1. Intended Use

The NADAL® RSV Test is a lateral flow chromatographic immunoassay for the qualitative detection of respiratory syncytial virus (RSV) antigens in human nasopharyngeal swab specimens. The test is intended for use as an aid in the diagnosis of RSV infections in patients showing typical symptoms of a respiratory tract infection (see section 12 'Limitations'). The test procedure is not automated and requires no special training or qualification. The NADAL® RSV Test is designed for professional use only.

2. Introduction and Clinical Significance

RSV is the most common respiratory virus in infants and young children. It infects virtually all children by the age of two. In most children, the virus causes symptoms resembling those of the common cold. In infants born prematurely and/or with chronic lung disease, RSV can cause a severe or even life-threatening disease. Prior to the introduction of Synagis®, RSV disease resulted annually in over 125,000 hospitalisations. There was a high mortality risk amongst approximately 2% of infants. Symptoms of RSV can include fever and a runny nose as well as other, more serious symptoms like coughing, difficult and rapid breathing or wheezing.

3. Test Principle

The NADAL® RSV Test enables the detection of RSV antigens through the visual interpretation of colour development on the internal test strip. Monoclonal antibodies specific to RSV antigens are immobilised in the test line region (T) of the membrane. During the test, the specimen reacts with the anti-RSV antibodies which are conjugated to coloured particles and precoat onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there are a sufficient number of RSV antigens in the specimen, a coloured line will develop in the test line region (T) of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result. The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® RSV test cassettes
 - 20 extraction tubes containing buffer (approximately 300 µL each)*, incl. dropper caps
 - Additional material provided according to 93/42/EEC:
Due to possible supply shortages of COVID-19 related accessory medical products, the swab manufacturer might change. Therefore, the supplied swabs are from one of the manufacturers listed below.
- a) 20 sterile swabs, CE 0197

Jiangsu Changfeng Medical Industry Co. Ltd.,
Touqiao Town, Guangling District, Yangzhou,
Jiangsu 225109 China (authorised EU
representative: Llins Service & Consulting GmbH,
Obere Seegasse 34/2, 69124 Heidelberg,
Germany)

- b) 20 sterile swabs, CE 0197

CITOTEST LABWARE MANUFACTURING Co. Ltd.,
No. 339 Beihai West Road, Haimen, 226100
Jiangsu, P.R. China
(authorised EU representative: WellKang Ltd,
Enterprise Hub, NW Business Complex,
1 Beraghmore Rd., Derry, BT48 8SE, Northern
Ireland)

- c) 20 sterile swabs, CE 0197

Hangzhou Yiguoren Biotechnology Co. Ltd., Room
402, Building 2, No. 2628, Yuhangtang Road,
Cangqian Street, Yuhang District, Hangzhou,
311121 Zhejiang, China (authorised EU
representative: Zoustech S.L., Paseo de la
Castellana 141, 28049 Madrid, Spain)

- d) 20 sterile swabs, CE 2797

Puritan Medical Products Company LLC,
31 School Street, Guilford, Maine 04443-0149
USA (authorised EU representative: EMERGO
EUROPE, The Hague, The Netherlands)

- 1 reagent holder
- 1 package insert

*containing the following preservative: ProClin™ 300:
<0.03%.

No hazard labelling for ProClin™ 300 is required according to Regulation (EC) № 1272/2008 CLP. Concentrations are below the exemption threshold of <0.03%.

5. Additional Materials Required

- Timer

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the entire instructions for use prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.

- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Do not use swabs from damaged pouches.
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

Nasopharyngeal specimen:


- It is important to obtain as much secretion as possible.
- Insert the swab into the nostril that presents the most secretion under visual inspection, parallel to the palate (not upwards), until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx.
- Gently rub and roll the swab. Leave the swab in place for several seconds to absorb secretions.
- Slowly remove the swab whilst rotating it. Specimens can be collected from both nostrils using the same swab, but it is not necessary to collect specimens from both sides if the tip is saturated with fluid from the first collection.

Note:


- Swab specimens should be tested immediately after collection. For best test performance, use freshly collected specimens.
- Do not use specimens that are obviously contaminated with blood, as it may interfere with the flow of specimens and lead to inaccurate test results.
- To ensure optimal performance, use the swabs supplied in the test kit.

9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Hold a specimen collection tube containing buffer and labelled with the patient or control identification vertically. Peel off the aluminum foil from the tube.
 
2. Insert the swab with the collected specimen into the tube. Swirl the swab and squeeze it at least 3 times by
 


compressing the wall of the extraction tube against the swab to extract the antigens contained in the swab.

3. Place the tube into the designated area of the reagent holder.
 
4. Let the solution stand for 1 minute.
 
5. Remove the swab, pressing it firmly against the wall of the tube to release as much liquid as possible. Dispose of the swab in accordance with guidelines for the handling of infectious agents.
6. Attach a dropper cap to the tube.
 
7. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.
 
8. Place the test cassette on a clean and level surface.
9. Invert the tube and add 3 drops (approximately 100 µl) of the extracted solution to the specimen well (S) of the test cassette.
 
10. Start the timer.
 
11. Wait for the coloured line(s) to appear. Read the test result after 15 minutes. Do not interpret the result after more than 20 minutes.

10. Result Interpretation

Positive:

A coloured line develops in the control line region (C) and another coloured line develops in the test line region (T).

Note: The colour intensity in the test line region (T) may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region (T) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in the control line region (C). No line develops in the test line region (T).

Invalid:

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

The internal procedural control is included in the test cassette: A coloured line appearing in the control line region (C) is considered an internal positive procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® RSV Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of RSV antigens in human nasopharyngeal swab specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of RSV antigens can be determined with this qualitative test.
- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').
- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods (e.g. cell culture or a molecular assay) is recommended. A negative result does not at any time preclude the possibility of an RSV infection.
- The NADAL® RSV Test only detects the presence of RSV antigens in specimens and should not be used as the sole criterion for a diagnosis of an RSV infection.
- Positive and negative predictive values are highly dependent on prevalence. The local prevalence should be taken into consideration when interpreting diagnostic test results.
- Positive results do not preclude co-infections with other pathogens.
- As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® RSV Test was evaluated using clinical nasopharyngeal swab specimens in comparison to PCR.

The results are presented in the following table:

NADAL® RSV Test		PCR		
		Positive	Negative	Total
	Positive	63	0	63
	Negative	6	133	139
	Total	69	133	202

Diagnostic sensitivity: 91.3% (82.3% - 96.0%)*

Diagnostic specificity: >99.9% (97.2% - 100.0%)*

Overall agreement: 97.0% (93.7% - 98.6%)*

*95% confidence interval

Analytical performance

Detection limit

The 95% detection limit of the NADAL® RSV Test for RSV subtype A strains is 93.6 TCID₅₀/mL and RSV subtype B strains is 11.7 TCID₅₀/mL, and was determined via the dilution of inactivated RSV.

Measuring range

No adverse effect on T-line formation (prozone effect) was observed when testing specimens containing an RSV antigen concentration as high as 4.68×10^4 TCID₅₀/mL for RSV subtype A strains and as high as 1.17×10^5 TCID₅₀/mL for RSV subtype B strains. Thus, the measuring range of the test is between at least 93.6 TCID₅₀/mL and 4.68×10^4 TCID₅₀/mL for RSV subtype A strains and between at least 11.7 TCID₅₀/mL and 1.17×10^5 TCID₅₀/mL for RSV subtype B strains.

Analytical specificity

Cross-reactivity study

RSV negative swabs were spiked with 10 µL of the following potentially cross-reacting substances and tested using the NADAL® RSV Test:

Substance	Concentration
<i>Streptococcus pneumoniae</i>	1×10^6 CFU/mL
Influenza A/California/04/2009	1×10^6 TCID ₅₀ /mL
Influenza B/Washington/02/2019	1×10^5 TCID ₅₀ /mL
Epstein-Barr virus	1×10^6 copies/mL
Group A <i>Streptococcus</i>	1×10^6 CFU/mL

No cross-reactivity with the specimens was observed when tested using the NADAL® RSV Test.

Interference study

RSV negative swabs, swabs spiked with 10 µL of RSV A positive control and swabs spiked with 10 µL of RSV B positive control were spiked with 10 µL of the following potentially interfering substances and showed no interference with the NADAL® RSV Test:

Substance	Concentration	Substance	Concentration
4-acetamidophenol	10 mg/mL	Blood, type A	2% v/v
Acetylsalicylic acid	20 mg/mL	Blood, type B	2% v/v
Albuterol	0.083 mg/mL	Blood, type AB	2% v/v
Amantadine	500 ng/mL	Blood, type O	2% v/v
Mometasone	500 ng/mL	Dexamethasone	10 mg/mL
Mouthwash Listerine®	5% v/v	Dextromethorphan	10 mg/mL
Mouthwash Lion	5% v/v	Diphenhydramine HCl	5 mg/mL
Ascorbic acid chewable tablets	5% w/w	Ibuprofen	10 mg/mL
Beclomethasone	500 ng/mL	Mucin	1 mg/mL

Substance	Concentration	Substance	Concentration
Nasal spray	5% v/v	Oseltamivir	500 ng/mL
Oxymetazoline	0.05 mg/mL	Phenylephrine	1 mg/mL
Zanamivir	1 mg/mL	Tobramycin	500 ng/mL
Benzocaine throat spray Cepacol®	5% v/v		

Precision

Repeatability

Repeatability was established by testing 20 replicates of negative, low, medium and high positive nasopharyngeal swab specimens using 3 lots of the NADAL® RSV tests. >99% of the specimens were correctly identified (20/20 correct tests per concentration, 95% confidence interval: 95.5% - 100%). The NADAL® RSV Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 5 replicates of negative, low, medium and high positive nasopharyngeal swab specimens. Testing was performed by 5 operators using 3 independent NADAL® RSV test lots at 3 different sites on 5 separate days. >99% of the specimens were correctly identified (375/375 correct tests per concentration, 95% confidence interval: 99% - 100%). The NADAL® RSV Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the case of any serious incidents related to the performance of the NADAL® RSV Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C.B.,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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c) 20 écouvillons stériles, CE 0197

Hangzhou Yiguoren Biotechnology Co. Ltd., Room
402, Building 2, No. 2628, Yuhangtang Road,
Cangqian Street, Yuhang District, Hangzhou,
311121 Zhejiang, China

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15. Bibliographie

- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
- Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
- Hall, C,B,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
- Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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- a) 20 hisopos estériles CE 0197

Jiangsu Changfeng Medical Industry Co. Ltd.,
Touqiao Town, Guangling District, Yangzhou,

Jiangsu 225109 China (representante autorizado de la UE: Llins Service & Consulting GmbH, Obere Seegasse 34/2, 69124 Heidelberg, Germany)

- b) 20 hisopos estériles CE 0197

CITOTEST LABWARE MANUFACTURING CO. Ltd.,
No. 339 Beihai West Road, Haimen, 226100
Jiangsu, P.R. China
(representante autorizado de la UE: WellKang Ltd, Enterprise Hub, NW Business Complex, 1 Beragham Rd., Derry, BT48 8SE, Northern Ireland)

- c) 20 hisopos estériles CE 0197

Hangzhou Yiguoren Biotechnology Co. Ltd., Room 402, Building 2, No. 2628, Yuhangtang Road, Cangqian Street, Yuhang District, Hangzhou, 311121 Zhejiang, China (representante autorizado de la UE: Zoustech S.L., Paseo de la Castellana 141, 28049 Madrid, Spain)

- d) 20 hisopos estériles CE 2797

Puritan Medical Products Company LLC.
31 School Street, Guilford, Maine 04443-0149
USA (representante autorizado UE: EMERGO EUROPE, The Hague, The Netherlands)

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- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.
- Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89 : 1443-1447.
- Hall, C.B.J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N.Engl. J. Med. 308 : 1443-1447.
- Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. J.Clin. Microbiol. 257 : 763-767.

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15. Referencias

- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66 : 291-300.

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a) 20 tamponi sterili, CE 0197

Jiangsu Changfeng Medical Industry Co. Ltd.,
Touqiao Town, Guangling District, Yangzhou,
Jiangsu 225109 China (authorised EU
representative: Llns Service & Consulting GmbH,

Obere Seegasse 34/2, 69124 Heidelberg,
Germany)

b) 20 tamponi sterili CE 0177

CITOTEST LABWARE MANUFACTURING CO. Ltd.,
No. 339 Beihai West Road, Haimen, 226100
Jiangsu, P.R. China
(authorised EU representative: WellKang Ltd,
Enterprise Hub, NW Business Complex, 1
Beraghmore Rd., Derry, BT48 8SE, Northern
Ireland)

c) 20 tamponi sterili CE 0197

Hangzhou Yiguo Biotechnology Co. Ltd., Room
402, Building 2, No. 2628, Yuhangtang Road,
Cangqian Street, Yuhang District, Hangzhou,
311121 Zhejiang, China (authorised EU
representative: Zoustech S.L., Paseo de la
Castellana 141, 28049 Madrid, Spain)

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15. Bibliografia

- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.

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(representante autorizado da UE: WellKang Ltd, Enterprise Hub, NW Business Complex, 1 Beraghmore Rd., Derry, BT48 8SE, Northern Ireland)

c) 20 zaragatoas estéreis, CE 0197

Hangzhou Yiguoren Biotechnology Co. Ltd., Room
402, Building 2, No. 2628, Yuhangtang Road,
Cangqian Street, Yuhang District, Hangzhou,
311121 Zhejiang, China (representante autorizado da UE: Zoustech S.L., Paseo de la Castellana 141, 28049 Madrid, Spain)

d) 20 zaragatoas estéreis, CE 2797

Puritan Medical Products Company LLC,
31 School Street, Guilford, Maine 04443-0149
USA (representante autorizado da UE: EMERGO EUROPE, The Hague, The Netherlands)

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15. Referências

- Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee corza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. Am. J. Hyg. 66 : 291-300.
- Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. J. Am. Med. Assoc. 176 : 647-653.
- Hall, C.B., R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. J. Pediatr. 89 : 1443-1447.
- Hall, C.B., J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. N.Engl. J. Med. 308 : 1443-1447.
- Taber, L.H.V., Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. Pediatrics 72 : 613-618.
- Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. J.Clin. Microbiol. 257 : 763-767.

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15. Reference

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
3. Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
4. Hall, C.B,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
5. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
6. Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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CITOTEST LABWARE MANUFACTURING Co. Ltd.,

No. 339 Beihai West Road, Haimen, 226100
Jiangsu, P.R. China
(valtuutettu EU-edustaja: WellKang Ltd,
Enterprise Hub, NW Business Complex, 1
Beraghmore Rd., Derry, BT48 8SE, Northern
Ireland)

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Hangzhou Yiguoren Biotechnology Co. Ltd., Room
402, Building 2, No. 2628, Yuhangtang Road,
Gangqian Street, Yuhang District, Hangzhou,
311121 Zhejiang, China (valtuutettu EU-edustaja:
Zoustech S.L., Paseo de la Castellana 141, 28049
Madrid, Spain)

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Puritan Medical Products Company LLC,
31 School Street, Guilford, Maine 04443-0149
USA (valtuutettu EU-edustaja: EMERGO EUROPE,
The Hague, The Netherlands)

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15. Lähteet

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
3. Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
4. Hall, C.B.,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
5. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
6. Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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15. Referenser

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
3. Hall, C.B, R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
4. Hall, C.B,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and WW.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
5. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
6. Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

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d) 20 sterile vattpinner, CE 0197

Puritan Medical Products Company LLC,
31 School Street, Guilford, Maine 04443-0149
USA (autorisert EU-representant: EMERGO
EUROPE, Haag, Nederland)

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15. Referanser

1. Chanock, R.M., and L. Findberg. 1957. Recovery from infants with respiratory illness of a virus related to chimpanzee coryza agent (CCA). II. Epidemiologic aspects of infection in infants and young children. *Am. J. Hyg.* 66 : 291-300.
2. Chanock, R.M., H.W. Kim, A.J. Vargosko, A. Deleva, K.M. Johnson, C.Cumming, and R.H. Parrott. 1961. Respiratory syncytial virus. I. Virus recovery and other observations during 1960 outbreak of bronchiolitis, pneumonia, and minor respiratory diseases in children. *J. Am. Med. Assoc.* 176 : 647-653.
3. Hall, C.B. R.G. Douglas, and J.M. Geiman. 1976. Respiratory syncytial virus infections in infants: quantitation and duration of shedding. *J. Pediatr.* 89 : 1443-1447.
4. Hall, C.B.,J.T. McBride, E.E. Walsh, D.M. Bell, C.L. Gala, S.Hildreth, L.G. Teneyck, and W.W.J. Hall. 1983. Aerosolized ribavirin treatment of infants with respiratory syncytial virus infection. *N.Engl. J. Med.* 308 : 1443-1447.
5. Taber, L.H.V, Knight, B.E. Gilbert, H.W. McClung, S.Z. Wilson, H.J. Norton, J.M. Thurson, W.H. Gordon, R.L. Atmar and W.R. Schlaudt. 1983. Ribavirin aerosol treatment of bronchiolitis associated with respiratory syncytial virus infection in infants. *Pediatrics* 72 : 613-618.
6. Ahluwalia, G.J. Embree, P. McNicol, B.Law, and G.W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunoabsorbent assay. *J.Clin. Microbiol.* 257 : 763-767.

Rev.0, 2022-05-25 Tr. ag.

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico- diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń
	Bevollmächtigter in der Europäischen Gemeinschaft	Authorised representative in the European Community	Mandataire dans l'Union Européenne	Representante autorizado en la Comunidad Europea	Mandatario nella Comunità Europea	Upoważniony przedstawiciel we Wspólnocie Europejskiej
	Bei beschädigter Verpackung nicht verwenden	Do not use if package is damaged	Ne pas utiliser si l'emballage est endommagé	No utilizar si el envase está dañado	Non utilizzare se la confezione è danneggiata	Nie używać w przypadku uszkodzonego opakowania

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøyte
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> -diagnostikkään tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur-begränsning	Temperatuurlimiet	Temperatur-begrænsning	Temperatur begrensnng
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov éčíslo	Luettelonumero	Listnummer	Catalogus nummer	Best il lingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for<n> tester
	Mandatário na Comunidade Europeia	Zplnomocněný zástupce v Evropském společenství	Euroopan yhteisössä toimiva valtuutettu edustaja	Auktoriserad representant inom Europeiska gemenskapen	Gemachtigde in de Europese Gemeenschap	Angiver den autoriserede repræsentant i Den Europæiske Union	Autorisert representant i EU
	Não utilizar se a embalagem estiver danificada	Nepoužívejte, je-li obal poškozen	Käyttö kielletty jos pakkaus on vaurioitunut	Får inte användas om förpackningen är skadad	Niet gebruiken indien de verpakking beschadigd is	Må ikke anvendes, hvis emballagen er beskadiget	Skal ikke brukes hvis pakningen er skadet

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0

Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0

Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29

Free Tel: 0800 291 565

Fax: +49 290 10-50

Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18

Free Tel –UK: 0808 234 1237

Free Tel – IRE: 1800 555 080

Fax: +49 290 10-50

France:

France Tel: 0800 915 240

France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720

Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82

Belgium Fax: 0800 747 07

Luxembourg

Lux, Tel: 800 211 16

Lux, Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759

Free Tel: 900 938 315

Fax: +49 941 290 10-50

Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34

Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44

Free Tel: 00 800 491 15 95

Fax: +49 941 290 10-50

Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735

Tel, Verde: 800 849 230

Fax: +49 941 290 10-50

Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600

Free Tel: 0800 0222 890

Fax: +31 70 30 30 775

Free Fax: 0800 024 9519

Nordic countries:

Denmark

Tel: +31 703075 605

Free Tel: 808 887 53

Finland

Tel: +31 703075 606

Free Tel: 0800 918 263

Free Fax: 0800 918 262

Norway

Tel: +31 703075 605

Free Tel: 800 16 731

Sweden

Tel: +31 703075 605

Free Tel: 020 79 09 06

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Tel: +49 2841 99820-0 • Fax: +49 2841 99820-1